

IH *Talent*

World's first versatile
adhesive applicator kit


1 click 1 dose


primequal
SURGICAL



IH Talent

World's first versatile adhesive applicator kit



Intended Use

Primequal™ Talent IH™ M.I.S. adhesive applicator kit is a groundbreaking cyanoacrylate microdosing delivery system developed to allow a safer, faster and solid bonding of surgical mesh.

Inspiring award-winning design

Talent IH™ is the first of a next-generation of adhesive applicator designed for surgeons to look into new, simple and reliable solutions to improve patient care and safety.

Versatile

Talent IH™ is the world's first versatile adhesive applicator kit. This is because it can be used with different glues*, it can fix multiple meshes and based on its cumulated skills, it will adapt to any situation.

No training required

Talent IH™ enables surgeons to perform a precise delivery immediately without any training, saving time and money and providing increased comfort and precision for themselves.

Safer

Talent IH™ innovative patented technology controls the precision and recurrence of the micro dose simply and comfortably. It has the unique ability to produce a **Contactless Delivery**, administering a dose safely and precisely without being in contact with the mesh or the tissue.

Faster

Talent IH™ automatic single button operation enables a faster procedure, including **Rapid Multipoint Delivery** under time pressure. Micro dosing delivers thin layers which is an accelerating factor for polymerization of the cyanoacrylates.

Automatic dose precision

AutomInject™ patented technology enables **Talent IH™** to deliver a precise micro dose each time the lever is pressed. It is a fully automatic **dose per dose** delivery from the first to the last dose. **AutomInject™** produces a « click » sound after every dose, enabling the surgeon to effectively monitor the correct dosage: **1 click = 1 dose**. The treatment is therefore greatly simplified to a click based procedure.

No air bubble

Talent IH™ has been designed to ensure that no air bubble can disturb the procedure.

No wastage

Talent IH™ should be filled with the required dose for the surgery. Liquid path has been optimised to minimise dead space and wastage.


Overpressure safety

Talent IH™ is equipped with an overpressure safety feature that will avoid administering glue incorrectly.

Compatibility

Talent IH™ is compatible with most cyanoacrylates*. Filling **Talent IH™** is identical to filling a standard syringe, simple and fast.

*1 : Refer to Primequal™ list of compatible glues
Talent IH™ is a single use and disposable delivery system. Standard and regulatory references: European Directive 93 / 42 / EEC, ISO Standard 13485 and ISO Standard 14971. The manufacturer reserves the right to modify the product without prior notice. Primequal™ and Talent™ are Trademarks of Primequal SA, all rights reserved.



reddot design award



100% Swiss Engineering
100% Swiss Design
100% Swiss Made
100% Swiss Quality



Safer Faster Injectors

Rue des Pierres du Niton 17 • 1207 Geneva • Switzerland
Tel + 41 22 354 05 50 • Fax + 41 22 354 05 51 • www.primequal.com • info@primequal.com

Pharmaceutical
Innovation of the year
Mention: «Ease of Use»

